

[KO Validated] B-Catenin Rabbit Mab For Wb, Ip, If-P, Ihc-P, Elisa - P35222

Item Number: CM0003509



Introduction

Rabbit monoclonal antibody against human β -Catenin (P35222), KO validated for WB, IP, IF-P, IHC-P, and ELISA. Reacts with human, mouse, rat. Ideal for Wnt signaling and cell adhesion studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] β -Catenin Rabbit mAb
Remarks/Alias	EVR7; CTNNB; MRD19; NEDSDV; armadillo; in
Species	Human
GeneID	1499
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human beta Catenin (P35222).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35222
Protein Weight	85kDa
Shipping	Ice bag