


KO Validated Xiap Rabbit Polyclonal Antibody For Wb, Elisa - P98170

Item Number: CM0030934



Introduction

KO validated XIAP rabbit polyclonal antibody for WB, ELISA. Detects human, mouse, rat XIAP. XIAP is an apoptosis inhibitor and E3 ubiquitin ligase regulating NF-κB signaling, immunity, and cancer. Ideal for apoptosis and cell signaling research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] XIAP Rabbit pAb
Remarks/Alias	API3; ILP1; MIHA; XLP2; BIRC4; IAP-3; hIAP3; hIAP-3; AP
Species	Human
Gene ID (Human)	331
Gene ID	331
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 198-497 of human XIAP (NP_001158.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P98170
Protein Weight	57kDa
Shipping	Ice bag