

[KO Validated] Ythdf2 Rabbit Mab - Q9Y5A9

Item Number: CM0006545



Introduction

Validated rabbit monoclonal antibody against YTHDF2, an m6A reader protein (62 kDa) regulating mRNA stability. Suitable for WB, IHC-P, IP, ELISA with cross-reactivity to human, mouse, rat.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] YTHDF2 Rabbit mAb
Remarks/Alias	DF2; CAHL; HGRG8; NY-REN-2; YTHDF2
Species	Human
Gene ID (Human)	51441
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y5A9
Protein Weight	62 kDa
Shipping	Ice bag