

Kpnb1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q14974

Item Number: CM0025296



Introduction

KPNB1 rabbit polyclonal antibody targets importin subunit beta-1, a key nuclear transport receptor. Validated for WB, IHC-P, IF/ICC, ELISA in human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	KPNB1 Rabbit pAb
Remarks/Alias	IMB1; IPO1; IPOB; Impnb; NTF97; KPNB1
Species	Human
Gene ID (human)	3837
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 300-400 of human KPNB1 (NP_002256.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14974
Protein Weight	97kDa
Shipping	Ice bag