

Kras Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P01116

Item Number: CM0018516



Introduction

High-quality KRAS rabbit polyclonal antibody for Western blot, immunofluorescence, and ELISA. Detects human, mouse, and rat KRAS. Essential tool for Ras-MAPK signaling and cancer research.

[Learn More](#)

Parameter	Value
Product Name	KRAS Rabbit pAb
Remarks/Alias	NS; NS3; OES; CFC2; RALD; K-Ras; KRAS1; KRAS2; RASK2; KI-RAS; C-K-RAS; K-RAS2A; K-RAS2B; K-RAS4A; K-RAS4B; K-Ras 2; 'C-K-RAS; c-Ki-ras; c-Ki-ras2; KRAS
Species	Human
GeneID (human)	3845
GeneID	3845
Immunogen	Synthetic peptide corresponding to amino acids 30-189 of human KRAS (NP_203524.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01116
Protein Weight	22kDa
Shipping	Ice bag