

Kras+Hras+Nras Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P01111 / P01112 / P01116

Item Number: CM0005168



Introduction

Pan-Ras (KRAS, HRAS, NRAS) Rabbit Monoclonal Antibody validated for WB, IF/ICC, ELISA in human, mouse, rat. Detects ~18kDa band, targeting the conserved N-terminus; ideal for Ras-MAPK signaling and cancer research.

[Learn More](#)

Parameter	Value
Product Name	KRAS+HRAS+NRAS Rabbit mAb
Remarks/Alias	NS6; CMNS; KRAS; NCMS; ALPS4; N-ras; NRAS1; KRAS+HRAS+NRAS
Species	Human
GeneID (Human)	4893 / 3265 / 3845 (NRAS / HRAS / KRAS)
Immunogen	Synthetic Peptide corresponding to aa 1-100 of human NRAS (P01111)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01111 / P01112 / P01116
Protein Weight	18kDa
Shipping	Ice bag