


L3MbtI3 Rabbit Monoclonal Antibody For Wb, Ip, Elisa - Q96Jm7

Item Number: CM0003299



Introduction

Rabbit monoclonal antibody against human L3MBTL3 (MBT-1), a chromatin regulator involved in Notch signaling and protein degradation. Validated for WB, IP and ELISA in human samples.

[Learn More](#)

Parameter	Value
Product Name	L3MBTL3 Rabbit mAb
Remarks/Alias	MBT1; MBT-1; L3MBTL3
Species	Human
Gene ID (Human)	84456
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-255 of human L3MBTL3 (NP_115814.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human
SWISS	Q96JM7
Protein Weight	88kDa
Shipping	Ice bag