

Lactic Acid-Histone H3-K18 Rabbit Polyclonal Antibody (Q16695/P68431) For Wb, Db, Elisa

Item Number: CM0000086



Introduction

Rabbit polyclonal antibody targeting lactylated lysine 18 of human Histone H3, suitable for WB, DB, and ELISA. Detects human, mouse, and rat samples with broad predicted reactivity. Ideal for epigenetic and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	Lactic acid-Histone H3-K18 Rabbit pAb
Remarks/Alias	H3/A; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FA; H3C10; H3C11; H3C12; HIST1H3A; Lactic acid-Histone H3-K18
Species	Human
GeneID	8290/8350
Immunogen	A synthetic lactylated peptide around K18 of human Histone H3 (NP_003520.1)
Source	Rabbit
Category	Other Modified Antibodies
Application	WB, DB, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag