


Lactylated Histone H3 (K14) Rabbit Polyclonal Antibody For Wb/Elisa (Q16695/P68431)

Item Number: CM0000089



Introduction

Rabbit polyclonal antibody specific for lactylated Histone H3 at lysine 14 (K14). Validated for Western blot (WB) and ELISA. Detects human, mouse, and rat Histone H3 lactylation, with wide range predicted cross-reactivity. Targeting SWISS Q16695/P68431.

[Learn More](#)

Parameter	Value
Product Name	Lactic acid-Histone H3-K14 Rabbit pAb
Remarks/Alias	H3/A; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FA; H3C10; H3C11; H3C12; HIST1H3A; Lactic acid-Histone H3-K14
Species	Human
GeneID	8290/8350
Immunogen	Synthetic peptide: A synthetic lactylated peptide around K14 of human Histone H3 (NP_003520.1).
Source	Rabbit
Category	Other Modified Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15kDa
Shipping	Ice bag