


Lactylated Histone H3-K18 Rabbit Monoclonal Antibody (Q16695/P68431) For Epigenetic Research

Item Number: CM0000080



Introduction

Rabbit monoclonal antibody against lactylated histone H3 at lysine 18 (SWISS: Q16695/P68431). Validated for WB and ELISA with broad cross-reactivity (human, mouse, rat). Suitable for studying histone lactylation and metabolic gene regulation.

[Learn More](#)

Parameter	Value
Product Name	Lactic acid-Histone H3-K18 Rabbit mAb
Remarks/Alias	H3/A; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FA; H3C10; H3C11; H3C12; HIST1H3A; Lactic acid-Histone H3-K18
Species	Human
GeneID	8290/8350
Immunogen	A synthetic lactylated peptide around K18 of human Lactic acid-Histone H3-K18 (NP_003520.1).
Source	Rabbit
Category	Other Modified Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15kDa
Shipping	Ice bag