


Lamin B1 Rabbit Mab For Wb, Ihc-P - P20700

Item Number: CM0002431



Introduction

Lamin B1 Rabbit monoclonal antibody for WB and IHC-P, validated in human, mouse, and rat. A key nuclear lamina component, ideal for chromatin organization and laminopathy research.

[Learn More](#)

Parameter	Value
Product Name	Lamin B1 Rabbit mAb
Remarks/Alias	LMN; ADLD; LMN2; LMNB; MCPH26
Species	Human
GeneID (Human)	4001
GeneID	4001
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 450-550 of human Lamin B1 (P20700).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P
Cross Reactivity	Human, Mouse, Rat
SWISS	P20700
Protein Weight	66kDa
Shipping	Ice bag