


Lamin B2 Rabbit Monoclonal Antibody For Wb - Q03252

Item Number: CM0008168



Introduction

Rabbit monoclonal antibody targeting Lamin B2 for Western blotting, IHC, IF, and ELISA. Cross-reacts with human, mouse, rat. Suitable for nuclear lamina and cell cycle studies.

[Learn More](#)

Parameter	Value
Product Name	Lamin B2 Rabbit mAb
Remarks/Alias	EPM9; LMN2; LAMB2; MCPH27; Lamin B2
Species	Human
GeneID (Human)	84823
GeneID	84823
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q03252
Protein Weight	70kDa
Shipping	Ice bag