


Lamtor5 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O43504

Item Number: CM0026745



Introduction

LAMTOR5 Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Targets Ragulator complex protein LAMTOR5 (HBXIP), involved in mTORC1 signaling and apoptosis.

[Learn More](#)

Parameter	Value
Product Name	LAMTOR5 Rabbit pAb
Remarks/Alias	XIP; HBXIP; LAMTOR5
Species	Human
GeneID (Human)	10542
GeneID	10542
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-173 of human LAMTOR5 (NP_006393.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43504
Protein Weight	10kDa
Shipping	Ice bag