


Lc3B Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Gzq8

Item Number: CM0017964



Introduction

LC3B Rabbit Polyclonal Antibody (pAb) validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat LC3B. Ideal for autophagy and mitochondrial quality control studies.

[Learn More](#)

Parameter	Description
Product Name	LC3B Rabbit pAb
Remarks/Alias	LC3B; ATG8F; MAP1LC3B-a; MAP1A/1BLC3
Species	Human
Gene ID	81631
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human LC3B (NP_073729.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9GZQ8
Protein Weight	15kDa
Shipping	Ice bag