

Ldl Receptor (Ldlr) Rabbit Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P01130

Item Number: CM0002034



Introduction

High-quality rabbit monoclonal antibody against human LDL Receptor (LDLR), validated for WB, IF/ICC, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for cholesterol metabolism and viral entry research.

[Learn More](#)

Parameter	Value
Product Name	LDL Receptor (LDLR) Rabbit mAb
Remarks/Alias	FH; FHC; FHCL1; LDLCQ2; LDL Receptor (LDLR)
Species	Human
GeneID (Human)	3949
GeneID	3949
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 761-860 of human LDL Receptor (LDLR) (NP_000518.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01130
Protein Weight	95kDa
Shipping	Ice bag