

Lect1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - 075829

Item Number: CM0026872

Introduction

Rabbit polyclonal antibody targeting human LECT1/Chondromodulin (CNMD). Validated for WB, IF/ICC, ELISA; reacts with human and mouse. Suitable for cartilage development, angiogenesis, and cardiac valve research.

[Learn More](#)

Parameter	Value
Product Name	LECT1 Rabbit pAb
Remarks/Alias	CHM1, CHM-I, LECT1, BRICD3, MYETS1
Species	Human
Gene ID (Human)	11061
Gene ID	11061
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 75-334 of human LECT1 (NP_008946.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	075829
Protein Weight	37kDa
Shipping	Ice bag