

Lgals2 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P05162

Item Number: CM0028143



Introduction

High-quality LGALS2 rabbit polyclonal antibody validated for Western blot, immunofluorescence, and ELISA in human and mouse. Recombinant immunogen, 15 kDa, UniProt P05162. Ideal for lectin research.

[Learn More](#)

Parameter	Specification
Product Name	LGALS2 Rabbit pAb
Remarks/Alias	HL14; LGALS2
Species	Human
GeneID (human)	3957
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-74 of human LGALS2 (NP_006489.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P05162
Protein Weight	15kDa
Shipping	Ice bag