

Lhx2 Rabbit Mab For Wb And Elisa - P50458

Item Number: CM0003429

Introduction

LHX2 Rabbit monoclonal antibody validated for Western blot and ELISA. Reacts with human, mouse, and rat LHX2, a transcriptional activator involved in lymphoid and neural development. Ideal for IVD research.

[Learn More](#)

Parameter	Value
Product Name	LHX2 Rabbit mAb
Remarks/Alias	LH2; hLhx2; LHX2
Species	Human
GeneID (Human)	9355
GeneID	9355
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human LHX2 (P50458).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50458
Protein Weight	44 kDa
Shipping	Ice bag