

## Limd1 Rabbit Pab - Q9Ugp4

Item Number: CM0019650

### Introduction

Rabbit polyclonal antibody against human LIMD1 for WB and ELISA. Reacts with human and mouse LIMD1 (72 kDa). Suitable for tumor suppression, cell signaling, and miRNA silencing studies.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | LIMD1 Rabbit pAb                                                                                                 |
| Remarks/Alias    | LIMD1                                                                                                            |
| Species          | Human                                                                                                            |
| GeneID (Human)   | 8994                                                                                                             |
| GeneID           | 8994                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human LIMD1 (NP_055055.1) |
| Source           | Rabbit                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                            |
| Application      | WB, ELISA                                                                                                        |
| Cross Reactivity | Human, Mouse                                                                                                     |
| SWISS            | Q9UGP4                                                                                                           |
| Protein Weight   | 72kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                          |