

Limk1 Rabbit Polyclonal Antibody For Wb And Elisa - P53667

Item Number: CM0019193



Introduction

High-quality LIMK1 rabbit polyclonal antibody for Western blot (WB) and ELISA. Reacts with human, mouse, and rat. Essential for actin filament dynamics and brain development.

[Learn More](#)

Parameter	Value
Product Name	LIMK1 Rabbit pAb
Remarks/Alias	LIMK; LIMK-1; LIMK1
Species	Human
Gene ID	3984
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 510-647 of LIMK1 (NP_002305.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P53667
Protein Weight	73kDa
Shipping	Ice bag