

Lin28A Rabbit Mab - Q9H9Z2

Item Number: CM0002184



Introduction

LIN28A Rabbit monoclonal antibody validated for WB, IF/ICC, ELISA. Targets human LIN28A, a key regulator of miRNA biogenesis and pluripotency. Ideal for stem cell and developmental research.

[Learn More](#)

Parameter	Value
Product Name	LIN28A Rabbit mAb
Remarks/Alias	CSDD1; LIN28; LIN-28; ZCCHC1; lin-28A; LIN28A
Species	Human
GeneID	79727
Immunogen	A synthetic peptide corresponding to amino acids 100-209 of human LIN28A (NP_078950.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q9H9Z2
Protein Weight	23 kDa
Shipping	Ice bag