

Lipin 1 Rabbit Mab For Wb And Elisa - Q14693

Item Number: CM0010097

Introduction

Rabbit monoclonal antibody against human Lipin 1, validated for WB and ELISA. Cross-reacts with mouse and rat. Suitable for lipid metabolism, adipocyte differentiation, and mitochondrial dynamics research.

[Learn More](#)

Parameter	Value
Product Name	Lipin 1 Rabbit mAb
Remarks/Alias	PAP1; Lipin 1
Species	Human
Gene ID (Human)	23175
Gene ID	23175
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 101-200 of human Lipin 1 (Q14693).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14693
Protein Weight	99kDa
Shipping	Ice bag