


Lipocalin-2/Ngal Rabbit Pab - P80188

Item Number: CM0023747



Introduction

High-quality rabbit polyclonal antibody targeting human Lipocalin-2/NGAL, validated for WB, IHC-P, IF-P, and ELISA. Cross-reacts with mouse. Ideal for studying iron metabolism, innate immunity, and renal diseases.

[Learn More](#)

Parameter	Value
Product Name	Lipocalin-2/NGAL Rabbit pAb
Remarks/Alias	p25, 24p3, MSFI, NGAL, Lipocalin-2/NGAL
Species	Human
Gene ID	3934
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-198 of human Lipocalin-2/NGAL (NP_005555.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P80188
Protein Weight	23 kDa
Shipping	Ice bag