

Llgl2 Rabbit Polyclonal Antibody For Western Blot And Elisa - Llgl2

Item Number: CM0014010

Introduction

LLGL2 Rabbit Polyclonal Antibody (pAb) validated for Western blot (WB) and ELISA detection of mouse and rat LLGL2 protein (~114kDa). Suitable for research applications.

[Learn More](#)

Parameter	Value
Product Name	LLGL2 Rabbit pAb
Remarks/Alias	Llglh2; 9130006H11Rik; LLGL2
Species	Mouse
Gene ID (Human)	3993
Gene ID	217325
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 938-1027 of mouse LLGL2 (NP_663413.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
Protein Weight	114kDa
Shipping	Ice bag