

Lox Rabbit Mab For Wb, Ihc-P, Ip, Elisa - P28300

Item Number: CM0003861

Introduction

LOX (Lysyl oxidase) rabbit monoclonal antibody validated for WB, IHC-P, IP, ELISA. Reacts with human, mouse, rat LOX. Ideal for extracellular matrix and fibrosis research.

[Learn More](#)

Parameter	Value
Product Name	LOX Rabbit mAb
Remarks/Alias	AAT10; LOX
Species	Human
Gene ID (Human)	4015
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 318-417 of human LOX (P28300)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28300
Protein Weight	47kDa
Shipping	Ice bag