

Lpin2 Rabbit Polyclonal Antibody For Western Blot And Elisa - Q92539

Item Number: CM0018656

Introduction

LPIN2 rabbit polyclonal antibody validated for Western blot and ELISA. Reacts with human and mouse LPIN2, a phosphatidate phosphatase involved in lipid metabolism. Suitable for metabolic and disease research.

[Learn More](#)

Parameter	Value
Product Name	LPIN2 Rabbit pAb
Remarks/Alias	LPIN2; lipin 2
Species	Human
Gene ID	9663
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 180-400 of human LPIN2 (NP_055461.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q92539
Protein Weight	99kDa
Shipping	Ice bag