


Lrp1 Rabbit Polyclonal Antibody For Wb, Elisa - Q07954

Item Number: CM0013533



Introduction

LRP1 Rabbit Polyclonal Antibody validated for WB and ELISA, reacting with human and mouse LRP1, a 505 kDa endocytic receptor critical in lipid metabolism and TAU spread. Supports cardiovascular and neuroscience research.

[Learn More](#)

Parameter	Details
Product Name	LRP1 Rabbit pAb
Remarks/Alias	APR; KPA; LRP; A2MR; CD91; APOER; LRP1A; TGFBRS; IGFBP3R; IGFBP-3R; IGFBP3R1; LRP1
Species	Human
Gene ID (Human)	4035
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-168 of LRP1 (NP_002323.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q07954
Protein Weight	505 kDa
Shipping	Ice bag