

Lsamp Rabbit Polyclonal Antibody For Wb, Elisa - Q13449

Item Number: CM0015513

Introduction

This LSAMP rabbit polyclonal antibody is suitable for Western blot and ELISA, detecting human, mouse, and rat LSAMP. Useful for investigating neuronal guidance and hippocampal mossy fiber development.

[Learn More](#)

Parameter	Value
Product Name	LSAMP Rabbit pAb
Remarks/Alias	LAMP; IGLON3; LSAMP
Species	Human
Gene ID (Human)	4045
Gene ID	4045
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 219-315 of human LSAMP (NP_002329.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13449
Protein Weight	37kDa
Shipping	Ice bag