

Lumican (Lum) Rabbit Pab For Wb, If/Icc, Elisa - P51884

Item Number: CM0027355

Introduction

Lumican (LUM) rabbit polyclonal antibody validated for WB, IF/ICC, ELISA. Specific to human, mouse, rat lumican. Recombinant immunogen. Suitable for extracellular matrix and corneal research.

[Learn More](#)

Parameter	Value
Product Name	Lumican (LUM) Rabbit pAb
Remarks/Alias	LDC; SLRR2D; Lumican (LUM)
Species	Human
GeneID (Human)	4060
GeneID	4060
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51884
Protein Weight	38kDa
Shipping	Ice bag