

Ly-6A/E (Sca-1) Rabbit Pab - P05533

Item Number: CM0022513

Introduction

High-quality rabbit polyclonal antibody targeting mouse Ly-6A/E (Sca-1), validated for Western blot and ELISA applications. Ideal for T-cell activation and stem cell research.

[Learn More](#)

Parameter	Value
Product Name	Ly-6A/E (Sca-1) Rabbit pAb
Remarks/Alias	TAP; Sca1; Sca-1; Ly-6A.2; Ly-6A/E; Ly-6E.1; Ly-6A/E (Sca-1)
Species	Mouse
Gene ID	110454
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of mouse Ly6a. (NP_034868.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P05533
Protein Weight	14kDa
Shipping	Ice bag