

Ly-6C Rabbit Pab For Wb, Elisa - P0Cw02

Item Number: CM0019330

Introduction

Ly-6C Rabbit Polyclonal Antibody validated for Western blot and ELISA, targeting mouse lymphocyte antigen 6C1 (GPI-anchored membrane protein). Suitable for immune cell research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	Ly-6C Rabbit pAb
Remarks/Alias	Ly6c; Ly-6C; Ly-6C1
Species	Mouse
GenelD	17067 (Mouse)
Immunogen	Recombinant protein of mouse Ly-6C
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P0CW02
Protein Weight	14KDa
Shipping	Ice bag