

Ly-6G Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P35461

Item Number: CM0007191

Introduction

High-performance Ly-6G rabbit monoclonal antibody for WB, IF-P, IHC-P, and ELISA. Validated in mouse, targeting the 14 kDa lymphocyte antigen 6G (P35461) for bone marrow and immunology research.

[Learn More](#)

Parameter	Value
Product Name	Ly-6G Rabbit mAb
Remarks/Alias	Gr1; Gr-1; Ly-6G
Species	Mouse
GeneID	546644
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-105 of mouse Ly-6G (NP_001297367.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse
SWISS	P35461
Protein Weight	14kDa
Shipping	Ice bag