


LYAR Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9NX58

Item Number: CM0021555



Introduction

LYAR Rabbit Polyclonal Antibody (Q9NX58) validated for WB, IHC-P, IF/ICC, ELISA. Reacts with Human, Mouse, Rat. Ideal for studying rRNA processing and transcription regulation roles of LYAR.

[Learn More](#)

Parameter	Value
Product Name	LYAR Rabbit pAb
Remarks/Alias	ZLYAR; ZC2HC2; LYAR
Species	Human
GeneID (human)	55646
GeneID	55646
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 10-160 of human LYAR (NP_060286.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NX58
Protein Weight	44kDa
Shipping	Ice bag