

Lysozyme (Lyz) Rabbit Polyclonal Antibody - P61626

Item Number: CM0016845



Introduction

Lysozyme (LYZ) rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human and mouse lysozyme C, a bacteriolytic enzyme. Ideal for studying antimicrobial activity and amyloidosis.

[Learn More](#)

Parameter	Value
Product Name	Lysozyme (LYZ) Rabbit pAb
Remarks/Alias	LZM; LYZF1; Lysozyme (LYZ)
Species	Human
GeneID (human)	4069
GeneID	4069
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Lysozyme (Lysozyme (LYZ)) (NP_000230.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P61626
Protein Weight	17kDa
Shipping	Ice bag