


Mad1/Mad1L1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y6D9

Item Number: CM0014408



Introduction

High-quality rabbit polyclonal antibody against human MAD1 (MAD1L1), validated for WB, IF/ICC, and ELISA. Ideal for spindle assembly checkpoint and cell cycle regulation research.

[Learn More](#)

Parameter	Value
Product Name	MAD1/MAD1L1 Rabbit pAb
Remarks/Alias	MAD1; MVA7; PIG9; TP53I9; TXBP181; MAD1/MAD1L1
Species	Human
GeneID (human)	8379
GeneID	8379
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 360-700 of human MAD1L1 (NP_001013859.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q9Y6D9
Protein Weight	83kDa
Shipping	Ice bag