


Mad2B/Mad2L2 Rabbit Mab For Wb, If/Icc, Elisa - Q9Ui95

Item Number: CM0011920



Introduction

High-quality MAD2B/MAD2L2 (REV7) rabbit monoclonal antibody for Western Blot, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Ideal for DNA repair and cell cycle studies.

[Learn More](#)

Parameter	Value
Product Name	MAD2B/MAD2L2 Rabbit mAb
Remarks/Alias	REV7; FANCV; MAD2B; POLZ2; MAD2B/MAD2L2
Species	Human
GeneID (Human)	10459
GeneID	10459
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 100-200 of human MAD2B/MAD2B/MAD2L2 (Q9Ui95).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Ui95
Protein Weight	24 kDa
Shipping	Ice bag