


Madcam1 Rabbit Pab For Wb, Ihc-P, Elisa - Q13477

Item Number: CM0027991



Introduction

MADCAM1 rabbit polyclonal antibody suitable for WB, IHC-P, ELISA. Reacts with human, mouse, rat MADCAM1 (~40 kDa). Ideal for studying mucosal immunity and lymphocyte trafficking.

[Learn More](#)

Parameter	Value
Product Name	MADCAM1 Rabbit pAb
Remarks/Alias	MACAM1; MADCAM1
Species	Human
Gene ID (Human)	8174
Gene ID	8174
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human MADCAM1 (NP_570116.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13477
Protein Weight	40kDa
Shipping	Ice bag