

Mafb Rabbit Pab For Wb And Elisa - Q9Y5Q3

Item Number: CM0017003

Introduction

MAFB Rabbit pAb is a polyclonal antibody targeting human transcription factor MafB (36kDa). Validated for WB and ELISA, cross-reacts with mouse and rat. Ideal for hematopoiesis and differentiation studies.

[Learn More](#)

Parameter	Value
Product Name	MAFB Rabbit pAb
Remarks/Alias	KRML; MCTO; DURS3; MAFB
Species	Human
GeneID (Human)	9935
GeneID	9935
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9Y5Q3
Protein Weight	36kDa
Shipping	Ice bag