

Magoh Rabbit Pab - P61326

Item Number: CM0030058



Introduction

MAGOH Rabbit Polyclonal Antibody for WB/ELISA applications. Recognizes human, mouse, rat MAGOH (17kDa). Ideal for research on pre-mRNA splicing, exon junction complex, and nonsense-mediated decay.

[Learn More](#)

Parameter	Value
Product Name	MAGOH Rabbit pAb
Remarks/Alias	MAGOH1; MAGOHA; MAGOH
Species	Human
GeneID (Human)	4116
GeneID	4116
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-146 of human MAGOH (NP_002361.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61326
Protein Weight	17kDa
Shipping	Ice bag