

Mammaglobin A Rabbit Mab - Q13296

Item Number: CM0002022

Introduction

High-quality Mammaglobin A Rabbit monoclonal antibody for Western Blot (WB), Immunohistochemistry (IHC-P), and ELISA applications. Validated with recombinant protein immunogen, human cross-reactivity. Ideal for breast cancer research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	Mammaglobin A Rabbit mAb
Remarks/Alias	MGB1; UGB2; PSBP1; Mammaglobin A
Species	Human
Gene ID (Human)	4250
Gene ID	4250
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-93 of human Mammaglobin A (Q13296).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q13296
Protein Weight	10kDa
Shipping	Ice bag