

Man1A2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O60476

Item Number: CM0016095

Introduction

High-quality MAN1A2 rabbit polyclonal antibody (UniProt O60476) validated for WB, IHC-P, and ELISA. Detects human and mouse mannosyl-oligosaccharide 1,2-alpha-mannosidase IB, a key enzyme in N-glycan trimming.

[Learn More](#)

Parameter	Value
Product Name	MAN1A2 Rabbit pAb
Remarks/Alias	MAN1B, MAN1A2
Species	Human
GeneID (Human)	10905
GeneID	10905
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 60-180 of human MAN1A2 (NP_006690.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	O60476
Protein Weight	73kDa
Shipping	Ice bag