

Map2 Rabbit Mab - P11137

Item Number: CM0005130

Introduction

Rabbit monoclonal MAP2 antibody validated for WB, IF, IHC, ELISA in human, mouse, rat. Targets MAP2 isoforms A, B, C, enabling neuronal cytoskeleton and dendrite morphology studies.

[Learn More](#)

Parameter	Value
Product Name	MAP2 Rabbit mAb
Remarks/Alias	MAP-2; MAP2A; MAP2B; MAP2C; MAP2
Species	Human
Gene ID (Human)	4133
Gene ID	4133
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-118 of human MAP2 (NP_002365.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-F, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11137
Protein Weight	200 kDa
Shipping	Ice bag