

Map7D1 Rabbit Pab For Wb, Elisa - Q3Kqu3

Item Number: CM0019437

Introduction

MAP7D1 Rabbit Polyclonal Antibody detects human and mouse MAP7 domain-containing protein 1. Validated for WB and ELISA, this antibody targets a microtubule-stabilizing protein involved in cell motility and neurite outgrowth.

[Learn More](#)

Parameter	Value
Product Name	MAP7D1 Rabbit pAb
Remarks/Alias	RPRC1; PARCC1; MAP7D1
Species	Human
GeneID (human)	55700
GeneID	55700
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-460 of human MAP7D1 (NP_060537.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q3KQU3
Protein Weight	93kDa
Shipping	Ice bag