


Mapre1 Rabbit Pab - Q15691

Item Number: CM0030779



Introduction

MAPRE1 Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA. Detects human, mouse, rat MAPRE1 (EB1) with recombinant immunogen. 30kDa protein, shipped on ice. Supports microtubule dynamics and cell division studies.

[Learn More](#)

Parameter	Value
Product Name	MAPRE1 Rabbit pAb
Remarks/Alias	EB1; MAPRE1
Species	Human
GeneID (Human)	22919
GeneID	22919
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 134-268 of human MAPRE1 (NP_036457.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15691
Protein Weight	30kDa
Shipping	Ice bag