

Marcksl1 Rabbit Polyclonal Antibody For Wb And Elisa - P49006

Item Number: CM0030080

Introduction

Rabbit polyclonal antibody targeting MARCKSL1 (MARCKS-like protein 1), validated for Western blot and ELISA in human, mouse, and rat samples. Useful for studying actin cytoskeleton regulation and cell migration.

[Learn More](#)

Parameter	Value
Product Name	MARCKSL1 Rabbit pAb
Remarks/Alias	F52, MLP, MRP, MLP1, MACMARCKS, MARCKSL1
Species	Human
Gene ID (Human)	65108
Gene ID	65108
Immunogen	Recombinant protein containing a sequence corresponding to amino acids 1-195 of human MARCKSL1 (NP_075385.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49006
Protein Weight	20kDa
Shipping	Ice bag