

Mast Cell Tryptase Rabbit Mab For Wb, Ihc-P, Ip, Elisa - Q15661

Item Number: CM0006055

Introduction

Rabbit monoclonal antibody against human mast cell tryptase (TPSAB1, Q15661), suitable for WB, IHC-P, IP, ELISA. Cross-reacts with human and mouse. Validated for research on mast cell degranulation and innate immunity.

[Learn More](#)

Parameter	Value
Product Name	Mast Cell Tryptase Rabbit mAb
Remarks/Alias	TPS1; TPS2; TPSB1; TPSB2; Tryptase-2
Species	Human
Gene ID (Human)	7177
Gene ID	7177
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15661
Protein Weight	31 kDa
Shipping	Ice bag