

Mat1A Rabbit Pab - Q00266

Item Number: CM0030057



Introduction

Rabbit polyclonal antibody targeting human MAT1A (44 kDa), validated for WB and ELISA, cross-reacts with mouse and rat. Suitable for studying S-adenosylmethionine biosynthesis in liver.

[Learn More](#)

Parameter	Value
Product Name	MAT1A Rabbit pAb
Remarks/Alias	MAT; SAMS; MATA1; SAMS1; MAT1A
Species	Human
GeneID (Human)	4143
GeneID	4143
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-395 of human MAT1A (NP_000420.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q00266
Protein Weight	44kDa
Shipping	Ice bag