

Mavs Rabbit Mab For Wb, Ihc-P, If/Icc, Elisa - Q7Z434

Item Number: CM0005454



Introduction

MAVS Rabbit Monoclonal Antibody validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human and mouse. Ideal for innate immunity research. MW 57kDa. - Q7Z434.

[Learn More](#)

Parameter	Value
Product Name	MAVS Rabbit mAb
Remarks/Alias	IPS1; VISA; IPS-1; CARDIF; MAVS
Species	Human
Gene ID (Human)	57506
Gene ID	57506
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 121-220 of human MAVS (NP_065797.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q7Z434
Protein Weight	57kDa
Shipping	Ice bag