


## Mavs Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q7Z434

Item Number: CM0027543



### Introduction

MAVS Rabbit Polyclonal Antibody (IPS1/VISA/CARDIF) validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human, mouse, and rat. Ideal for innate immunity and antiviral signal transduction studies.

[Learn More](#)

| Parameter        | Value                                                                                                                   |
|------------------|-------------------------------------------------------------------------------------------------------------------------|
| Product Name     | MAVS Rabbit pAb                                                                                                         |
| Remarks/Alias    | IPS1; VISA; IPS-1; CARDIF; MAVS                                                                                         |
| Species          | Human                                                                                                                   |
| GeneID (Human)   | 57506                                                                                                                   |
| GeneID           | 57506                                                                                                                   |
| Immunogen        | Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human MAVS (NP_065797.2). |
| Source           | Rabbit                                                                                                                  |
| Category         | Polyclonal Antibodies                                                                                                   |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                                                |
| Cross Reactivity | Human, Mouse, Rat                                                                                                       |
| SWISS            | Q7Z434                                                                                                                  |
| Protein Weight   | 57kDa                                                                                                                   |
| Shipping         | Ice bag                                                                                                                 |