


Mbd2 Rabbit Pab - Q9Ubb5

Item Number: CM0020761



Introduction

MBD2 Rabbit pAb is a polyclonal antibody raised in rabbit against human MBD2, suitable for WB and ELISA with cross-reactivity to mouse and rat. Ideal for chromatin remodeling and gene silencing research.

[Learn More](#)

Parameter	Value
Product Name	MBD2 Rabbit pAb
Remarks/Alias	DMTase; NY-CO-41; MBD2
Species	Human
GeneID (human)	8932
GeneID	8932
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 202-411 of human MBD2 (NP_003918.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9UBB5
Protein Weight	43kDa
Shipping	Ice bag